


Continuous Flow Extraction and Analysis of a Single Use Bioprocess System

Kevin Rowland, M.S. and Mark Jordi, Ph.D.
Analytical Services Division; 200 Gilbert St. Mansfield, MA 02048

Objective

Single use Bio-processing systems are attractive due to their efficiency and the potential to eliminate complex cleaning procedures and validations. This study was designed to demonstrate an analytical methodology tailored for effective determination of the extractable compounds from a model bioprocess system. A continuous flow extraction was performed in order to investigate the compounds extracted under the conditions most closely representing a use condition. Once compounds had been found and filtered, the effect of standard selection on the number of compounds falling above the analytical evaluation threshold (AET) was investigated.

Analytical Methodology



Model Bioprocess System

Bioprocess Bag

Multilayer Film
LDPE Liquid Contact
EVA Air Barrier Film

Filter

Nylon Filter
Polypropylene Housing
Ethylene propylene O-ring
Hydrophilic disk and holder

Tubing

Pharmaceutical grade
Thermoplastic Elastomer
Ultra-pure fluid Processing applications


Extraction Vehicle:

Saline- 137mM NaCl, 2.7mM KCl
10 mM Phosphate Buffer (pH 7.4),
300 mL




Extraction Conditions:
Continuous recirculation with a peristaltic pump; 72 hours; 37°C


Scouting Analyses



Elemental Analysis



Organic Identification



Analytical Evaluation Threshold

$$AET = \left(\frac{SCT}{\text{doses/day}} \right) \times \left(\frac{\text{labeled doses}}{\text{Container}} \right) = \frac{0.15 \mu\text{g}}{\text{day}} \times \frac{1 \text{ day}}{2 \text{ doses}} \times \frac{500 \text{ doses}}{\text{System}} = 37.5 \mu\text{g/Device}$$

Accounting for the extraction performed, the instrumental limit is: $37.5 \mu\text{g/Device} \div 300 \text{ mL} = 0.125 \mu\text{g/mL}$

Norwood D. The Analytical Evaluation Threshold (AET) and Its Relationship to Safety Thresholds, Leachables and Extractables Handbook (2012); 59-78. PQRI.org.


Number of compounds observed Above AET depends on Standard Selected



Analytical

Standard

Approx. AET Response Cut-off



Conclusions

A wide array of analytical techniques are required for thorough identification of system extractables. Extraction methodology effects the quantity of extractables observed. The analytical standard used to relate AET concentration to instrument response has a strong effect on the response cut-off. PQRI recommendation of adjusting the AET by the greater of 50% or the %RSD of the relative response factor database represents good compromise.